

SEC23 protein (STJ504891)

STJ504891

GENERAL INFORMATION

Product Type	Proteins
Short Description	SEC23-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10484
Gene Symbol	SEC23A
Uniprot ID	SC23A_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 700-765 on human SEC23A protein.
Immunogen Region	700-765
Specificity	
Immunogen Sequence	