

SEC24D protein (STJ504895)

STJ504895

GENERAL INFORMATION

Product Type	Proteins
Short Description	SEC24D-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9871
Gene Symbol	SEC24D
Uniprot ID	SC24D_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 860-920 on human SEC24D protein.
Immunogen Region	860-920
Specificity	
Immunogen Sequence	