

SEPT1 protein (STJ504896)
STJ504896**GENERAL INFORMATION**

Product Type	Proteins
Short Description	SEPT1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1731
Gene Symbol	SEPTIN1
Uniprot ID	SEPT1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-350 on SEPT1 protein.
Immunogen Region	300-350
Specificity	
Immunogen Sequence	