

Septin 2 protein (STJ504897)

STJ504897

GENERAL INFORMATION

Product Type	Proteins
Short Description	Septin 2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4735
Gene Symbol	SEPTIN2
Uniprot ID	SEPT2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 270-320 on human Septin-2 protein.
Immunogen Region	270-320
Specificity	
Immunogen Sequence	