

SEPT5 protein (STJ504900)
STJ504900**GENERAL INFORMATION**

Product Type	Proteins
Short Description	SEPT5-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5413
Gene Symbol	SEPTIN5
Uniprot ID	SEPT5_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-350 on SEPT5 protein.
Immunogen Region	300-350
Specificity	
Immunogen Sequence	