

MSF protein (STJ504904)
STJ504904**GENERAL INFORMATION**

Product Type	Proteins
Short Description	MSF-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10801
Gene Symbol	SEPTIN9
Uniprot ID	SEPT9_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-60 on SEPT9 protein.
Immunogen Region	1-60
Specificity	
Immunogen Sequence	