

SERHL protein (STJ504906)

STJ504906

GENERAL INFORMATION

Product Type	Proteins
Short Description	SERHL-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	253190
Gene Symbol	SERHL2
Uniprot ID	SEHL2_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on human serine hydrolase 1 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	