

SESN2 protein (STJ504909)

STJ504909

GENERAL INFORMATION

Product Type	Proteins
Short Description	SESN2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	83667
Gene Symbol	SESN2
Uniprot ID	SESN2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-80 on SESN2 protein.
Immunogen Region	1-80
Specificity	
Immunogen Sequence	