

SFRP2 protein (STJ504911)
STJ504911**GENERAL INFORMATION**

Product Type	Proteins
Short Description	SFRP2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6423
Gene Symbol	SFRP2
Uniprot ID	SFRP2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 225-295 on human SFRP2.
Immunogen Region	225-295
Specificity	
Immunogen Sequence	