

Soluble Guanylyl Cyclase Alpha 2 protein (STJ504914)

STJ504914

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Proteins |
| Short Description | Soluble Guanylyl Cyclase Alpha 2-protein if for use in research applications. |
| Applications | Immunodepletion, Immunocompetition |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | |
| Formulation | |
| Isotype | |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 66012 |
| Gene Symbol | Gucy1a2 |
| Uniprot ID | GCYA2_RAT |
| Immunogen | Synthetic peptide taken within amino acid region 680-730 on rat guanylate cyclase 1, soluble, alpha 2 protein. |
| Immunogen Region | 680-730 |
| Specificity | |
| Immunogen Sequence | |