

Soluble Guanylyl Cyclase beta 1 protein (STJ504915)

STJ504915

GENERAL INFORMATION

Product Type	Proteins
Short Description	Soluble Guanylyl Cyclase beta 1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25202
Gene Symbol	Gucy1b1
Uniprot ID	GCYB1_RAT
Immunogen	Synthetic peptide taken within amino acid region 175-225 on rat guanylate cyclase soluble subunit beta-1 protein.
Immunogen Region	175-225
Specificity	
Immunogen Sequence	