

SGK1 protein (STJ504916)

STJ504916

GENERAL INFORMATION

Product Type	Proteins
Short Description	SGK1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6446
Gene Symbol	SGK1
Uniprot ID	SGK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 40-90 on human SGK1 protein.
Immunogen Region	40-90
Specificity	
Immunogen Sequence	