

SGK3 protein (STJ504918)

STJ504918

GENERAL INFORMATION

Product Type	Proteins
Short Description	SGK3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100533105/23678
Gene Symbol	SGK3
Uniprot ID	SGK3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 440-490 on human Sgk3 protein.
Immunogen Region	440-490
Specificity	
Immunogen Sequence	