

SGLT2 protein (STJ504921)

STJ504921

GENERAL INFORMATION

Product Type	Proteins
Short Description	SGLT2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6524
Gene Symbol	SLC5A2
Uniprot ID	SC5A2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 580-630 on human SGLT2 protein.
Immunogen Region	580-630
Specificity	
Immunogen Sequence	