

SIGLEC10 protein (STJ504931)

STJ504931

GENERAL INFORMATION

Product Type	Proteins
Short Description	SIGLEC10-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	89790
Gene Symbol	SIGLEC10
Uniprot ID	SIG10_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 650-697 on human SIGLEC10 protein.
Immunogen Region	650-697
Specificity	
Immunogen Sequence	