

SIRT3 protein (STJ504934)
STJ504934**GENERAL INFORMATION**

Product Type	Proteins
Short Description	SIRT3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide: RDLVGPVLSPPRKDWQL corresponding to amino acids 203-221 of rat SIRT3 protein.
Immunogen Region	203-221
Specificity	
Immunogen Sequence	