

SIRT5 protein (STJ504935)
STJ504935**GENERAL INFORMATION**

Product Type	Proteins
Short Description	SIRT5-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	306840
Gene Symbol	Sirt5
Uniprot ID	SIR5_RAT
Immunogen	Synthetic peptide taken within amino acid region 260-310 on rat Sirtuin 5 protein.
Immunogen Region	260-310
Specificity	
Immunogen Sequence	