

Slingshot homolog 1 protein (STJ504938)

STJ504938

GENERAL INFORMATION

Product Type	Proteins
Short Description	Slingshot homolog 1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54434
Gene Symbol	SSH1
Uniprot ID	SSH1_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on Slingshot homolog 1 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	