

EAAT2 protein (STJ504939)

STJ504939

GENERAL INFORMATION

Product Type	Proteins
Short Description	EAAT2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6506
Gene Symbol	SLC1A2
Uniprot ID	EAA2_HUMAN
Immunogen	Synthetic peptide corresponding to C terminal amino acids 550-574 of Human EAAT2. Sequence: KVTLAANGKSADCSVEEEPWKREK
Immunogen Region	550-574
Specificity	
Immunogen Sequence	