

Smad4 protein (STJ504945)

STJ504945

GENERAL INFORMATION

Product Type	Proteins
Short Description	Smad4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4089
Gene Symbol	SMAD4
Uniprot ID	SMAD4_HUMAN
Immunogen	Synthetic peptides taken within amino acid region 400-450 and 502-552 on human SMAD4 protein.
Immunogen Region	400-450 and 502-552
Specificity	
Immunogen Sequence	