

SorLA protein (STJ504950)

STJ504950

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Proteins |
| Short Description | SorLA-protein if for use in research applications. |
| Applications | Immunodepletion, Immunocompetition |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | |
| Formulation | |
| Isotype | |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 20660 |
| Gene Symbol | Sor11 |
| Uniprot ID | SORL_MOUSE |
| Immunogen | Synthetic peptide taken within amino acid region 2100-2200 on human SORLA receptor protein. |
| Immunogen Region | 2100-2200 |
| Specificity | |
| Immunogen Sequence | |