

SOX10 protein (STJ504952)

STJ504952

GENERAL INFORMATION

Product Type	Proteins
Short Description	SOX10-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6663
Gene Symbol	SOX10
Uniprot ID	SOX10_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-450 on human Transcription Factor SOX-10 protein.
Immunogen Region	400-450
Specificity	
Immunogen Sequence	