

MASPIN protein (STJ504955)

STJ504955

GENERAL INFORMATION

Product Type	Proteins
Short Description	MASPIN-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5268
Gene Symbol	SERPINB5
Uniprot ID	SPB5_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 100-150 on human MASPIN protein.
Immunogen Region	100-150
Specificity	
Immunogen Sequence	