

**SRPX2 protein (STJ504963)**  
STJ504963**GENERAL INFORMATION**

<b>Product Type</b>	Proteins
<b>Short Description</b>	SRPX2-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

**PRODUCT PROPERTIES**

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">27286</a>
<b>Gene Symbol</b>	<a href="#">SRPX2</a>
<b>Uniprot ID</b>	<a href="#">SRPX2_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 430-470 on human SRPX2 protein.
<b>Immunogen Region</b>	430-470
<b>Specificity</b>	
<b>Immunogen Sequence</b>	