

STIM2 protein (STJ504964)

STJ504964

GENERAL INFORMATION

Product Type	Proteins
Short Description	STIM2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57620
Gene Symbol	STIM2
Uniprot ID	STIM2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 600-650 on human STIM2 protein.
Immunogen Region	600-650
Specificity	
Immunogen Sequence	