

TMEM73 protein (STJ504965)

STJ504965

GENERAL INFORMATION

Product Type	Proteins
Short Description	TMEM73-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution	
Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide corresponding to amino acid region 332-342. The selected peptide was post-synthetically modified to achieve highest antigenicity before coupling to a carrier protein using hetero bifunctional cross linker. Sequence: SFLLSSQEVLRHLLR
Immunogen	332-342
Region	
Specificity	
Immunogen	
Sequence	