

Stra6 protein (STJ504967)

STJ504967

GENERAL INFORMATION

Product Type	Proteins
Short Description	Stra6-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	20897
Gene Symbol	Stra6
Uniprot ID	STRA6_MOUSE
Immunogen	Synthetic peptide taken within third large intra cytoplasmic loop within amino acid region 450-500 on human Stra6 protein.
Immunogen Region	450-500
Specificity	
Immunogen Sequence	