

LYK5 protein (STJ504968)
STJ504968**GENERAL INFORMATION**

Product Type	Proteins
Short Description	LYK5-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	92335
Gene Symbol	STRADA
Uniprot ID	STRAA_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 50-100 on human STE-20 related kinase adapter protein alpha.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	