

SUR1 protein (STJ504970)
STJ504970**GENERAL INFORMATION**

Product Type	Proteins
Short Description	SUR1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25559
Gene Symbol	Abcc8
Uniprot ID	ABCC8_RAT
Immunogen	Synthetic peptide corresponding to positions 1560-1582.
Immunogen Region	1560-1582
Specificity	
Immunogen Sequence	