

SYDE1 protein (STJ504972)

STJ504972

GENERAL INFORMATION

Product Type	Proteins
Short Description	SYDE1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	85360
Gene Symbol	SYDE1
Uniprot ID	SYDE1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 680-735 on human SYDE1 protein.
Immunogen Region	680-735
Specificity	
Immunogen Sequence	