

Syndecan-1 protein (STJ504977)

STJ504977

GENERAL INFORMATION

Product Type	Proteins
Short Description	Syndecan-1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6382
Gene Symbol	SDC1
Uniprot ID	SDC1_HUMAN
Immunogen	C-epitope synthetic peptide taken within amino acid region 250-350 on human Syndecan-1 protein.
Immunogen Region	250-350
Specificity	
Immunogen Sequence	