

Synaptotagmin 1 protein (STJ504978)

STJ504978

GENERAL INFORMATION

Product Type	Proteins
Short Description	Synaptotagmin 1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6857
Gene Symbol	SYT1
Uniprot ID	SYT1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-60 on SYT1 protein.
Immunogen Region	1-60
Specificity	
Immunogen Sequence	