

NAK protein (STJ504981)

STJ504981

GENERAL INFORMATION

Product Type	Proteins
Short Description	NAK-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29110
Gene Symbol	TBK1
Uniprot ID	TBK1_HUMAN
Immunogen	Synthetic peptide taken within second intraceclular loop of human TBK1 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	