

ETS1 associated protein II protein (STJ504983)

STJ504983

GENERAL INFORMATION

Product Type	Proteins
Short Description	ETS1 associated protein II-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51567
Gene Symbol	TDP2
Uniprot ID	TYDP2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 330-360 on human Tyrosyl-DNA phosphodiesterase protein.
Immunogen Region	330-360
Specificity	
Immunogen Sequence	