

TDRD1 protein (STJ504984)

STJ504984

GENERAL INFORMATION

Product Type	Proteins
Short Description	TDRD1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	56165
Gene Symbol	TDRD1
Uniprot ID	TDRD1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 570-605 on TDRD1 protein.
Immunogen Region	570-605
Specificity	
Immunogen Sequence	