

TGM1 protein (STJ504990)

STJ504990

GENERAL INFORMATION

Product Type	Proteins
Short Description	TGM1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	74145
Gene Symbol	F13a1
Uniprot ID	F13A_MOUSE
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on Coagulation factor XIIIa chain precursor
Immunogen Region	
Specificity	
Immunogen Sequence	