

TIA1 protein (STJ504993)

STJ504993

GENERAL INFORMATION

Product Type	Proteins
Short Description	TIA1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: IHC: 1:50 ICC: 1:50 IF: 1:50 CFM: 1:50
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7072
Gene Symbol	TIA1
Uniprot ID	TIA1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 414-464 on human TIA1 protein.
Immunogen Region	414-464
Specificity	
Immunogen Sequence	