

TIPIN protein (STJ504996)
STJ504996**GENERAL INFORMATION**

Product Type	Proteins
Short Description	TIPIN-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54962
Gene Symbol	TIPIN
Uniprot ID	TIPIN_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 240-270 on mouse TIPIN protein.
Immunogen Region	240-270
Specificity	
Immunogen Sequence	