

TLR1 protein (STJ505001)

STJ505001

GENERAL INFORMATION

Product Type	Proteins
Short Description	TLR1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	21897
Gene Symbol	Tlr1
Uniprot ID	TLR1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 200-250 on mouse toll-like receptor 1 protein.
Immunogen Region	200-250
Specificity	
Immunogen Sequence	