

TLR11 protein (STJ505002)

STJ505002

GENERAL INFORMATION

Product Type	Proteins
Short Description	TLR11-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	239081
Gene Symbol	Tlr11
Uniprot ID	TLR11_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 430-460 on human TLR11 protein.
Immunogen Region	430-460
Specificity	
Immunogen Sequence	