

TLR4 protein (STJ505005)

STJ505005

GENERAL INFORMATION

Product Type	Proteins
Short Description	TLR4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7099
Gene Symbol	TLR4
Uniprot ID	TLR4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 970-1000 on human TLR4 protein.
Immunogen Region	970-1000
Specificity	
Immunogen Sequence	