

TLR5 protein (STJ505006)

STJ505006

GENERAL INFORMATION

Product Type	Proteins
Short Description	TLR5-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7100
Gene Symbol	TLR5
Uniprot ID	TLR5_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 195-225 on human TLR5 protein.
Immunogen Region	195-225
Specificity	
Immunogen Sequence	