

TLR13 protein (STJ505011)

STJ505011

GENERAL INFORMATION

Product Type	Proteins
Short Description	TLR13-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	279572
Gene Symbol	Tlr13
Uniprot ID	TLR13_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 400-480 of mouse TLR13 protein.
Immunogen Region	400-480
Specificity	
Immunogen Sequence	