

TMEM16A protein (STJ505016)

STJ505016

GENERAL INFORMATION

Product Type	Proteins
Short Description	TMEM16A-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55107
Gene Symbol	ANO1
Uniprot ID	ANO1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 900-950 on human Transmembrane protein 16A.
Immunogen Region	900-950
Specificity	
Immunogen Sequence	