

Timeless protein (STJ505018)

STJ505018

GENERAL INFORMATION

Product Type	Proteins
Short Description	Timeless-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	21853
Gene Symbol	Timeless
Uniprot ID	TIM_MOUSE
Immunogen	Synthetic peptide: SGTPRVHRKRFQIEDED corresponding to amino acids 1178-1197 on Mouse Timeless protein.
Immunogen Region	1178-1197
Specificity	
Immunogen Sequence	