

TOMM6 protein (STJ505021)

STJ505021

GENERAL INFORMATION

Product Type	Proteins
Short Description	TOMM6-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100188893
Gene Symbol	TOMM6
Uniprot ID	TOM6_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 24-74 on human Mitochondrial import receptor subunit TOM6 homolog.
Immunogen Region	24-74
Specificity	
Immunogen Sequence	