

Topoisomerase II alpha protein (STJ505023)

STJ505023

GENERAL INFORMATION

Product Type	Proteins
Short Description	Topoisomerase II alpha-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7153
Gene Symbol	TOP2A
Uniprot ID	TOP2A_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1220-1270 on human DNA topoisomerase (ATP-hydrolysing) protein.
Immunogen Region	1220-1270
Specificity	
Immunogen Sequence	