

TPCN1 protein (STJ505025)

STJ505025

GENERAL INFORMATION

Product Type	Proteins
Short Description	TPCN1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	53373
Gene Symbol	TPCN1
Uniprot ID	TPC1_HUMAN
Immunogen	Synthetic peptide taken within TMD7 and TMD8 on human TPC1 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	