

TPCN2 protein (STJ505026)

STJ505026

GENERAL INFORMATION

Product Type	Proteins
Short Description	TPCN2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	233979
Gene Symbol	Tpcn2
Uniprot ID	TPC2_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 680-731 on mouse TPC2 protein.
Immunogen Region	680-731
Specificity	
Immunogen Sequence	