

TPCN2 protein (STJ505027)

STJ505027

GENERAL INFORMATION

Product Type	Proteins
Short Description	TPCN2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	219931
Gene Symbol	TPCN2
Uniprot ID	TPC2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 205-255 on human TPC2 protein.
Immunogen Region	205-255
Specificity	
Immunogen Sequence	